

# A MALDI-TOF MS database for fast identification of *Vibrio* spp. potentially pathogenic in marine molluscs

VIVALDI Final conference  
Brest, 26-29 November

M. MOUSSA-POULY, M. A. TRAVERS, E. CAUVIN, B. THUILLIER, A. LE PIOUFFLE,  
O. LUCAS, M. TREILLES, C. GARCIA



# Our Objective...



Face to:

- Noticeable increase in bacterial diseases reports in marine molluscs  
↳ Bacteria of *Vibrio* genus : the genus the most involved during mortality events
- Difficulties in rapid identification of *Vibrio* spp.




Important to increase reactivity for monitoring program



**MALDI-TOF MS development**



**Quicker identification and characterization of *Vibrio* spp.**




# The MALDI-TOF MS (Matrix-assisted Laser Desorption/Ionisation Time of flight)

Mass spectrometry  
technique :


Bacterial proteins  
analysis



Detection




Protein spectrum :  
specific to a bacterium



Comparison  
Bacterial spectrum vs  
Spectra in a database

Identification



Acceleration

# How to develop MALDI-TOF MS ?

VibrioBase: A MALDI-TOF MS database for fast identification of *Vibrio* spp. that are potentially pathogenic in humans

René Erler<sup>a,\*</sup>, Antje Wichels<sup>a</sup>, Ernst-August Heinemeyer<sup>b</sup>, Gerhard Hauk<sup>c</sup>, Martin Hippelein<sup>d</sup>, Nadja Torres Reyes<sup>e</sup>, Gunnar Gerds<sup>a</sup>

<sup>a</sup> Alfred Wegener Institute for Polar and Marine Research, Biologische Anstalt Helgoland, Kurpromenade 201, 27498 Helgoland, Germany

<sup>b</sup> Governmental Institute of Public Health of Lower Saxony, Lüchtenburger Weg 24, 26603 Aurich, Germany

<sup>c</sup> Regional Office for Health and Social Affairs of Mecklenburg-Western Pomerania, Gertraudenstraße 11, 18057 Rostock, Germany

<sup>d</sup> University Medical Center Schleswig-Holstein, Central Facility: Medical Investigation Office and Hygiene, Brunswicker Straße 4, 24105 Kiel, Germany

<sup>e</sup> GFZ German Research Centre for Geosciences, Helmholtz Centre Potsdam, Telegrafenberg, 14473 Potsdam, Germany

Recent development of Maldi databases in marine sciences (Biotyper Bruker, Erler and al. 2015).

But these existing databases contain very few marine bacteria



Creation of a **Vibriobase**



A MALDI-TOF MS database  
Specific to Marine *Vibrio*

# How to generate this database ?


**PERTINENT  
BACTERIAL  
STRAINS  
SELECTION**



**BACTERIA  
STRAINS  
CULTURE**



**RAW  
SPECTRA  
ACQUISITION**



**SPECTRA  
SELECTION**

**REFERENCE  
SPECTRA  
CREATION**

**MALDI-TOF  
DATABASE**

# How to generate this database ?

SELECT  
PERTINENT  
BACTERIAL  
STRAINS

- Strains **from mortality** episodes (since 1988) in France
  - Isolated from oysters, mussels, clams, abalones, cockles
- = **121 strains**



*V. aestuarianus*



*V. harveyi*

*V. jasicida*

*V. rotiferianus*

*V. coralliilyticus*

*V. neptunius*

*V. ostreicida*


*V. bivalvicia*

*V. europaeus*

*V. tubiashii*

*V. orientalis*

*V. tapetis*




Other strains in *Harveyi*, *Mediterranei* and  
*Orientalis* clades

# How to generate this database ?

**SELECT  
PERTINENT  
BACTERIAL  
STRAINS**

**BACTERIA  
STRAINS  
CULTURE**

- Culture on three different media :
  - Zobell
  - 1,5% salt TSA
  - Marine Agar
- By three collaborating labs (Labeo, Laboceia and Qualyse)
- To increase the diversity of conditions to obtain strong results




# How to generate this database ?

**SELECT  
PERTINENT  
BACTERIAL  
STRAINS**

**BACTERIA  
STRAINS  
CULTURE**

**RAW SPECTRA  
ACQUISITION**


- Generation of 24 raw spectra / strain / medium
- 72 spectra / strain



# How to generate this database ?

## SPECTRA SELECTION

- Use of different methods of spectra analysis :
  - Visual method (Flexanalysis software)
  - Statistical methods (PCA clustering and CCI matrix)




- To keep the most specific spectra

# How to generate this database ?

SPECTRA  
SELECTION

REFERENCE  
SPECTRA  
CREATION

- **Reference spectrum** = MSP (main projection spectrum)
- **1MSP / strain** = compilation of selected spectra obtained in the three different media



# How to generate this database ?

## SPECTRA SELECTION


- All the MSP created constitute the Vibriodatabase
- This database contains 121 MSP

## REFERENCE SPECTRA CREATION

## MALDIT-TOF DATABASE

Vibrio aestuarianus francensis 12-128	1T3	64
Vibrio aestuarianus francensis 12-146	3T1	57
Vibrio aestuarianus francensis 12-147	1T1	69
Vibrio aestuarianus huître plate 04-091	1T1	64
Vibrio aestuarianus huître plate 04-091	2T1	58
Vibrio aestuarianus moule 16-055	1T1	63
Vibrio aestuarianus moule 16-055	5T1	69
Vibrio alginolyticus LMG 4409		67
Vibrio bivalvicia 14-105	2T1	60
Vibrio bivalvicia 14-105	3T1	69
Vibrio bivalvicia 14-106	1T1	62
Vibrio bivalvicia 14-106	49T1	66
Vibrio bivalvicia CECT 8855		61
Vibrio bivalvicia CECT 8856		62
Vibrio brasiliensis DSMZ 17184		60
Vibrio campbelli CIP 70.67		64
Vibrio corallilyticus 06-209		45
Vibrio corallilyticus 06-210		67
Vibrio corallilyticus 06-211		64
Vibrio corallilyticus 06-221		70
Vibrio corallilyticus 06-225		66
Vibrio corallilyticus 09-123	1T5	47
Vibrio corallilyticus LMG 19607		53
Vibrio corteqadensis DSM 102805		70
Vibrio europaeus 07-038	2T2	68
Vibrio europaeus 07-108	T1	60
Vibrio europaeus 07-110	T1	70
Vibrio europaeus 07-112	T1	54
Vibrio europaeus 07-115	T2	47
Vibrio europaeus 07-116	T1	68
Vibrio europaeus 07-117	T1	66
Vibrio europaeus 07-119	T1	69
Vibrio europaeus 07-120	T1	52
Vibrio europaeus 07-121	1T1	67
Vibrio europaeus CECT 8136		65
Vibrio harveyi 08-075	3T2	61
Vibrio harveyi 08-077	4T2	68
Vibrio harveyi carchariae 04-018	1T1	66
Vibrio harveyi carchariae 04-092	2T4	69
Vibrio harveyi carchariae 04-107	1T1	64


# Step 1 : test of the created MSP with existing databases



MSP of species present in databases



96%  
of MSP well-identified




MSP of species not present in databases



No identification or miss-  
identification for all strains

## Step2 : Create a dendogramm with our newly created MSP



# Step 3: MALDI-TOF Vibriobase validation



## BLIND TEST OF 100 VIBRIO STRAINS

## RAW SPECTRA ACQUISITION

## DATABASE TEST


- ✓ Analysed by three French collaborators (three different media)
- ✓ Two protocols were tested **Direct transfer** (bacteria colony) and **protein extract**



- ✓ 3 raw spectra obtained / strain : **direct transfer**
- ✓ 3 raw spectra obtained / strain: **protein extract**

- ✓ Obtained spectra were tested with :
  - The newly created Vibriobase
- ✓ To be validated

Logarithmic score value




# Validation : blind test with 100 new bacterial strains

## Lab 1

### Direct transfer

94 %



### Protein extract

94 %


- With direct transfer and protein extract : majority of tested strains well-identified
- The mismatching occurred mainly with :
  - V. europaeus* / *V. tubiashii*
  - V. hepatarius* / *V. orientalis*
  - V. owensii* / *V. jasicida*

# Validation : blind test with 100 new bacterial strains

Lab 2

Direct transfer


Protein extract



- With direct transfer and protein extract : majority of tested strains well-identified
- The mismatching occurred mainly with :
  - V. europaeus* / *V. tubiashii*
  - V. hepatarius* / *V. orientalis*
  - V. owensii* / *V. jasicida*

# Validation : blind test with 100 new bacterial strains


## Lab 3



- With direct transfer and protein extract : majority of tested strains well-identified
- The mismatching occurred mainly with :
  - V. europaeus* / *V. tubiashii*
  - V. hepatarius* / *V. orientalis*
  - V. owensii* / *V. jasicida*

# Focus on identification limits...

- High matching score (>2,3) of :
  - *V. europaeus* with *V. tubiashii*
  - *V. hepatarius* with *V. orientalis*
  - *V. owensii* with *V. jasicida*
- Species very close regarding MALDI-TOF biotyper Tree
- *V. europaeus* and *V. tubiashii* recently dissociated and were classified as the same species
- Just one strain of *V. owensii* in the database




Dendrogram on the basis of the MSP spectra (Maldi Biotyper)

# Conclusion

- The MALDI-TOF database created contains 121 strains of *Vibrio* specific to marine molluscs
- A majority of *Vibrio* species were well-identified with a good specificity
- MALDI-TOF Database can be used as routine diagnostic for rapid identification
- The direct transfer method can be used for greater spread

MSP Name
<i>Vibrio bivalvicia</i> 14-106 49T1 QUALYSE 22
<i>Vibrio bivalvicia</i> 14-106 49T1 56
<i>Vibrio bivalvicia</i> 14-106 49T1 66
<i>Vibrio bivalvicia</i> CECT 8855
<i>Vibrio bivalvicia</i> CECT 8855 LABEO
<i>Vibrio bivalvicia</i> CECT 8855 LABEO 20
<i>Vibrio bivalvicia</i> CECT 8855 LABOCEA
<i>Vibrio bivalvicia</i> CECT 8855 LABOCEA 20
<i>Vibrio bivalvicia</i> CECT 8855 LABOCEA 23
<i>Vibrio bivalvicia</i> CECT 8855 QUALYSE
<i>Vibrio bivalvicia</i> CECT 8855 61
<i>Vibrio bivalvicia</i> CECT 8855 66
<i>Vibrio bivalvicia</i> CECT 8856
<i>Vibrio bivalvicia</i> CECT 8856 LABEO
<i>Vibrio bivalvicia</i> CECT 8856 LABEO 22
<i>Vibrio bivalvicia</i> CECT 8856 LABOCEA
<i>Vibrio bivalvicia</i> CECT 8856 LABOCEA 20
<i>Vibrio bivalvicia</i> CECT 8856 QUALYSE
<i>Vibrio bivalvicia</i> CECT 8856 QUALYSE 21
<i>Vibrio bivalvicia</i> CECT 8856 QUALYSE 23
<i>Vibrio bivalvicia</i> CECT 8856 62
<i>Vibrio bivalvicia</i> CECT 8856 64
<i>Vibrio coralliilyticus</i> 06-209
<i>Vibrio coralliilyticus</i> 06-209 LABEO
<i>Vibrio coralliilyticus</i> 06-209 LABEO 21
<i>Vibrio coralliilyticus</i> 06-209 LABOCEA
<i>Vibrio coralliilyticus</i> 06-209 LABOCEA 20
<i>Vibrio coralliilyticus</i> 06-209 QUALYSE
<i>Vibrio coralliilyticus</i> 06-209 45
<i>Vibrio coralliilyticus</i> 06-209 65
<i>Vibrio coralliilyticus</i> 06-209 68
<i>Vibrio coralliilyticus</i> 06-210
<i>Vibrio coralliilyticus</i> 06-210 LABEO
<i>Vibrio coralliilyticus</i> 06-210 LABOCEA
<i>Vibrio coralliilyticus</i> 06-210 QUALYSE 23



Real progress compared to time-consuming classical methods  
to identify correctly *Vibrio*

# Perspectives

- To make this database freely accessible
- To extend this database by including bacteria of the *Splendidus* clade
- To characterize *Vibrio* species at population level : identification of specific peaks spectrum to discriminate vir + / -

The final objective :

## THE VIBRIODATABASE



- ✓ Sequencing and MALDI-TOF databases for characterization
- ✓ To recommend relevant method for species identification (housekeeping genes, genome sequencing or MALDI-TOF)
- ✓ To be able to provide well-characterized strains to other European labs in order to identify quicker new genotype in Europe

LNR-LRUE



# Thank you for your attention !

---

## Acknowledgments :

LEMAR-Plouzané – Jean-Louis Nicolas (*V. coralliiolyticus*, *V. aestu.*)

LEMAR-Plouzané – Christine Paillard (*V. tapetis*)

PFOM-Roscoff – Frédérique Le Roux (*Splendidus* and *Harveyi* clades)

IHPE-Montpellier – Delphine Destoumieux (*Harveyi* clade)